

CERTIFIED®

A handy and portable spotting kit that can be used on any size spotting job.

SPOT & STAIN

Rescue Kit



Kit Contains

- Quart size bottles with pour tops
- Fresh Spot (W-1)
- Coffee / Tea (W-3)
- Digester
- Dr. Spot (Aerosol)
- White Towel
- Tamping Brush
- Spotting Chart

C361-004

SPOTTING CHART

Problem	Fresh Spot W-1	Coffee-Tea W-3	Digester	Dr-Spot
Adhesives				•
Alcohol Drinks		•		
Animal Glues			•	
Asphalt				•
Auto Grease				•
Beer		•		
Blood				•
Beverage Stains				•
Fresh Beverage	•			
Candy	•			
Chewing Gum				•
Clay/Dirt	•			
Chocolate		•		
Coffee/Tea		•		
Coffee/Tea w/cream	•	•		
Cola			•	
Cosmetics				•
Crayon				•
Cream	•		•	
Feces				•
Fruit Juices	•	•		
Greasy Foods				•
Ice Cream	•		•	
Inks (washable)	•			
Liquor	•			
Lipstick				•
Mildew			•	
Milk				•
Motor Oils				•
Mustard			•	
Paints (latex)			•	
Paint				•
Saliva		•	•	
Perfumes		•		
Shoe Polish				•
Soft Drinks		•		
Sugar		•		
Syrup		•		
Urine		•	•	
Tar				•
Vomit				•

All Certified's spotting agents are ready-to-use. All of the above problems can be removed with one or more spotters.

WHEN USING A SPOTTING KIT

- 1- Identify your spot, and reference on the chart above.
- 2- Pre-test products in inconspicuous location.
- 3- Apply spotter by blotting and tamping. Do not rub!
- 4- Always work from edge to center of spot.
- 5- Repeat on area until problem is gone.
- 6- Rinse with water. Do not over wet.

